



# Falls Prevention – Equipment installation & safe use

# **BEDSIDE COMMODE CHAIR**

These instructions are to be used in conjunction with the information that your Health Professional has discussed with you.





### WHAT IS A COMMODE CHAIR?


• A commode chair allows people who have mobility problems to toilet without having to get to the bathroom. This is particularly useful at night.

### **TO INSTALL**

- Most commodes have adjustable legs. The height should be adjusted to allow the user to sit and stand with ease. Ask your Occupational Therapist about the safest height.
- To adjust the height/leg length, push in the silver button to slide the internal leg piece up and down to the appropriate notch.

### HOW TO USE THE COMMODE CHAIR

- Raise the lid/seat of the commode.
- To sit on the seat move backwards until you can feel the commode behind your legs. Place your hands on the armrests for support then lower yourself slowly down.
- When getting up lean forward and push up on the arm rests









#### **PRECAUTIONS & SAFETY**

- Ensure that the commode chair is stable on the floor and cannot slide or tip over
- Ensure legs of the commode chair are adjusted to the correct height
- Ensure the toilet paper is situated within easy reach
- Ensure the pan is emptied regularly
- Ensure there is adequate space around the commode chair to enable safe movement on, off & around the commode

#### **MAINTENANCE & CLEANING**

- Ensure that you regularly clean the commode chair and bucket with a mild disinfectant.
- For a standard commode, lift the toilet seat and pick up the bucket using the handle.
- For a swing arm commode, your Occupational Therapist will show you how to swing the metal arm to the side, then pick up the bucket using the handle.
- Should the equipment appear damaged or corroded –do not use. It may require replacement.

Acknowledgement to: NSLHD and CCLHD Community Falls Prevention Committees. For further information scan this with your smart phone

Email:falls@cec.health.nsw.gov.auWeb:www.cec.health.nsw.gov.au



Clinical Excellence Commission©2016 Version 1 SHPN: 160355